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Cadastro ANVISA: N° 10216040034

Dental Molding Material
Alginate for Dental Impression
Type I Fast Setting

Indication for use:

The product is indicated for the making of dental arch impressions. Indicated in moldings for fixed and removable prostheses, for orthodontics and study models.

Warnings and Precautions:

- The inhalation of the powder may be dangerous to health.
- It can cause nausea and vomiting if swallowed.
- Exclusive use by the dental professional.
- Single-use product.

Instructions for Use:

- Use the following accessories: Measuring spoon for powder and graduated measuring cup with 3 water levels. One measure of powder (1 spoon = 10.5 g) and one measure of potable water (1 level = 18 ml) are enough to mold a hemiarch.
- Spatulate the powder and water until a homogenous mixture is obtained.
- Place the tray with ALGITEC PLUS in its proper position in the mouth, compressing it slightly.
- Make the muscular corrections and keep the tray still for 1 minute after the material starts to set.
- After setting, remove and rinse in running water.
- Work/spatulation time at 25°C: 35 seconds (± 05 sec.).
- Setting time: 1:30 min. (± 20 sec.)
- In case of delayed pouring and disinfection in the laboratory, after washing the impression, slightly shake it over the sink to remove excess water and place it in the humidity chamber with relative humidity of air at 90-100% or in a plastic bag without holes and without the addition of extra sources of humidity such as cotton, gauze or wet paper towels, place only the wet impression in the plastic bag and tie it or use zip-lock bags. This way the molding can wait for pouring for up to 100 hours.

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Packaging: Sachet of 20.5g, 410g and 454g.

Properties:

- Dustless;
- The formula was idealized to reduce the contraction by water loss (syneresis); the mold can also be stored for 100 hours in hermetically closed packaging before pouring, without suffering dimensional changes.
- Highly elastic;
- Complies with ISO 21563;
- Dosing according to the powder and water measurements;
- Chromatic - obtains a pink color during spatulation and returns to the white color after full gelation, which indicates the time of its removal from the mouth.

Composition:

Diatomite, Polyethylene Glycol, Calcium Sulfate, Potassium Alginate, Silicon Dioxide, Titanium Dioxide, Potassium Fluortitanate, Sodium Saccharin, Magnesium Oxide, Chlorhexidine Acetate, Tetrasodium Pyrophosphate and Artificial Aroma.

Storage: Keep away from heat and moisture. Open the packaging and transfer the content to a hermetic container. Always close the container immediately after use.

Handling: The product should be handled at temperature of approximately 20°C - 25°C at humidity below 70%.

Expiry: 2 years after date of manufacture.

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