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Cadastro ANVISA: Nº 10216040034

Dental Molding Material
Alginate for Dental Impression
Type II Normal Setting

Indication for use: The product is indicated for the making of dental arch impressions. Indicated in moldings for fixed and removable prostheses, for orthodontics and study models.

Warnings and Precautions:

- The inhalation of the powder may be dangerous to health.
- It can cause nausea and vomiting if swallowed.
- Exclusive use by the dental professional.
- Single-use product.

Instructions for Use:

- Use the following accessories: Measuring spoon for powder and graduated measuring cup with 3 water levels. One measure of powder (1 spoon = 10.5 g) and one measure of potable water (1 level = 18 ml) are enough to mold a hemiarch.
- Spatulate the powder and water until a homogenous mixture is obtained.
- Place the tray with DENCRIGEL in its proper position in the mouth, compressing it slightly.
- Make the muscular corrections and keep the tray still for 1 minute after the material starts to set.
- After setting, remove and rinse in running water.
- To obtain the best results, it is recommended to pour the plaster immediately after molding.
- Work/spatulation time at 25°C: 1:10 min. (± 10 sec.).
- Setting time: 2:00 min. (± 20 sec.)

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Packaging: Sachet of 20.5g, sachet of 410g and sachet of 454g.

Properties:

- Dustless
- Highly elastic
- Complies with ISO 21563
- Dosing according to the powder and water measurements
- Chromatic - obtains a pink color during spatulation and returns to the white color after full gelation, which indicates the time of its removal from the mouth.

Composition:

Diatomite, Polyethylene glycol, Calcium Sulfate, Potassium Alginate, Silicon Dioxide, Titanium Dioxide, Potassium Fluortitanate, Sodium Saccharin, Magnesium Oxide, Chlorhexidine Acetate, Tetrasodium Pyrophosphate and Artificial Aroma and Dye.

Storage: Keep away from heat and moisture. Open the packaging and transfer the content to a hermetic container. Always close the container immediately after use.

Handling: The product should be handled at temperature of approximately 20°C - 25°C at humidity below 70%.

Expiry: 2 years after date of manufacture.

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